

Colposcope YDJ II series

Product Description

ALLTION is proud to offer YDJ series of colposcope. Precise structure, new type of LED illumination system, adjustable balance spring arm and humanized design offer you excellent parfocal, large field view, precise and clear three dimensional images and ergonomic, fatigue-free operation.



YDJ series of colposcopes are the precise optical instruments specializing in gynecology examination and diagnose. Vagina lesion could be observed clearly which could not be seen by naked eyes through the instruments because of its magnification. It improved the precision of the diagnose. Clear image, wide field, even illumination, and convenient usage are the characteristics of YDJ series of colposcope. It is the necessary instruments for gynecology examination.

Features/Specifications

- 45° inclined binocular
- 16X Wide field eyepiece
- Smooth turret 5-step magnification pod
- Working distance 300MM
- Magnification 2.7X-30X
- Field diameter: 11mm-84mm
- Green and blue filter selected by a lever

**** Optional magnification/different objective combination provided**

- **The top 10W LED illumination**
- **20,000 hours LED lifespan**
- **30,000 LUX object illumination**

**** LED is lower energy consumption, longer lifespan, safer, more dependable than normal halogen illumination.**

- **Power supply voltage: 90V/60HZ-250V/50HZ**
- **DC12V output on the instruments and supply the power for CCD camera**

Optional accessories:

- ***20X Eyepiece***
- ***F=230mm objective***
- ***F=280mm objective***
- ***SLR digital camera adapter***
- ***Colposcope image collecting and processing system satisfy high-end observation and photograph. It is the best choice for digital camera.***