

# ***A New Star Is Born!***

*custo guard holter with custo wing –  
The best from four decades of experience  
in Holter ECG diagnostics!!*



## *custo guard holter with custo wing*

### *Modern Holter with novel patch electrodes*



custo guard holter is the new Holter ECG recorder that combines the benefits of both long-standing experience and continuous development. It excels in high patient comfort, is compact and lightweight, and comes without any disturbing cables or belts. With three channels, motion and pacemaker detection as well as an almost unlimited recording time of several months, the custo guard holter is the perfect solution for an excellent long-term ECG.

custo guard holter: The new all-in-one Holter monitor!

**Easy.**  
**Wireless.**  
**Reliable.**





## **Our benefits:**


- 3 channels with two independent leads
- Completely wireless and perfect for hygienic application
- Optimum recording quality with the new custo wing electrode
- Fast and stable device communication and data transfer
- Short battery charge times
- Pacemaker detection
- Position and motion detection
- Check for correct application through PC monitoring
- Bluetooth Low Energy for mobile applications
- LED status display for charge status and HR detection
- Protection class IP65: perfect for disinfection
- Very long recording time (several months)
- Running time up to 120 hours with only one battery
- Very small and lightweight with only 30 g
- Docking station custo base with one or six spots for device communication and as charging station

## **Various electrode applications are possible:**

- Special disposable electrode custo wing: Artefact reduction and signal optimisation due to decentral connections and fully covered adhesive surface – material made of breathable nonwoven fabric
- ECG cable guard 3 and ECG cable guard 4
- custo belt 3

## **Intelligent interface (guard holter viewer) in custo diagnostic**

- Display for recording stand-by mode and particular device states
- Various configurations and recording modes are possible



## Innovative system solutions for your diagnostics


Since 1982, our brand has been representing innovative system solutions in cardiopulmonary diagnostics. In cooperation with committed and competent sales and service partners, we offer professional solutions taking account of your individual organisational circumstances and IT structures. No matter if your installation is based on a stationary or mobile system, a single-user workstation or multi-site network solution, in a doctor's office or in a hospital – your modular custo med system will be a profitable investment for the future.

*Our modular platform custo diagnostic is the intelligent link between clinical engineering and administrative software.*



Visit us  
on YouTube!

[www.customed.de](http://www.customed.de)



custo med – your local dealer

Contact:  
custo med GmbH  
Maria-Merian-Str. 6  
85521 Ottobrunn  
Germany

Phone: +49 89 71098 00  
Fax: +49 89 71098 10  
Web: [www.customed.de](http://www.customed.de)  
E-mail: [info@customed.de](mailto:info@customed.de)

 **custo·med**  
EXCELLENCE IN DIAGNOSTICS