



**custo watch**  
**The new Holter Monitor**



**custo watch:**

**The result of consistent development work.**

The Holter ECG is recorded with a lightweight chest belt, transmitted wirelessly to the ECG watch and is saved there. It comes without disposable electrodes or cables but with a high degree of patient comfort and an attractive design.

The new custo watch provides noticeable and measurable benefits for both the patient and the doctor.

Reliable

Powerful

Ergonomic

Waterproof

Robust

Precise

**Chronometer**  
date, time

**Status**  
battery display, connection quality

**Info**  
ECG view, heartrate



**custo watch:**

**The advantages of a sophisticated system.**



- custo belt 3**
- Without adhesive electrodes
  - Novel belt technology for 3-channel ECG recording
  - Reusable
  - Machine washable
  - High wearing comfort
  - With certificate of biocompatibility



- custo guard 3**
- 3-channel ECG transmitter
  - Waterproof
  - Extremely lightweight
  - Energy-saving wireless transmission of ECG signal to custo watch



- custo watch**
- Signal testing of ECG channels
  - Adjustable patient display
  - Interchangeable bracelets
  - Waterproof
  - Extremely lightweight
  - Recording time: 72 h



- custo docking station**
- Charging of custo guard and custo watch
  - Fast data transfer of Holter ECG data to the PC
  - Compatible with Holter ECG evaluation software custo tera





## custo holter ABPM:

# Synchronous recording and display of Holter ECG and ABPM


### ABPM recorder custo screen 400

Beside the functions of custo screen 300, it additionally offers:

- Continuous saving of 3-channel ECG data from custo guard 3
- SD memory card for fast data download
- ABPM recording also possible without Holter ECG

### Software custo tera + custo screen

- Combined display of Holter ECG and ABPM



### Holter ECG transmitter custo guard 3

- 3-channel ECG recording
- Waterproof, shatterproof, extremely lightweight
- Energy-saving wireless transmission of ECG signal to the ABPM recorder custo screen 400

### 3-channel electrode belt custo belt 3

- Holter ECG without adhesive electrodes
- Novel belt technology for 3-channel ECG recording
- Reusable, machine washable
- High wearing comfort
- With certificate of biocompatibility



#### Contact:

custo med GmbH  
 Maria-Merian-Straße 6  
 85521 Ottobrunn  
 Germany

Phone: +49 89 71098 00  
 Fax: +49 89 71098 10  
 internet: [www.customed.de](http://www.customed.de)  
 e-mail: [info@customed.de](mailto:info@customed.de)

custo med - your local partner