



SIEMENS



Modern Healthcare, Meet Modern Technology.


ACUSON Freestyle Ultrasound System

Musculoskeletal


Meeting the evolving demands of healthcare, the ACUSON Freestyle™ ultrasound system empowers you to work the way you want using the most advanced wireless ultrasound technologies for ultrasound-guided procedures in musculoskeletal imaging.

For more information, visit
www.siemens.com/freestyle


Answers for life.



Botulinum toxin injection in lower leg for exertional compartment syndrome



Suprapatellar recess effusion



Transverse median nerve at the level of the wrist



Next Level Efficiency

Experience the freedom of wireless transducers with the ACUSON Freestyle system. Independence from cable management and decreased sterile sheath preparation time allow you to start your procedures faster. The system's imaging parameters can be managed remotely with integrated transducer controls for single operator use. Speed up procedures with automated uniform focusing, autofrequency, and auto TGC.



Value-based Care

Delivering quality and safety, the ACUSON Freestyle system provides the care your patients deserve with value-based technologies that improve patient outcomes and satisfaction. Pixelformer™ image processing architecture generates uniform imaging and excellent detail resolution to visualize the target and surrounding anatomy to effectively perform deep or superficial procedures. Enhanced needle visualization comes standard to quickly identify the needle both in- and out-of-plane. Wireless transducers maintain the sterile field for patient safety, and are fully submersible, including batteries, for high-level disinfection.



Innovation That Matters

Siemens provides the latest innovation to overcome clinical and operational challenges with the ACUSON Freestyle system. Cable-free, real-time scanning is made possible with proprietary ultra-wideband radio for rapid wireless data transfer. Bluetooth® radio powers the integrated transducer controls and a transducer alarm system. The ACUSON Freestyle system leads the way in healthcare with innovation that matters.

- Transducers are powered by rechargeable batteries for 90 minutes of continuous scanning.
- Enterprise-level wireless networking comes standard so you can send patient study lists to PACS or network storage.
- Ability to view and share images and clips from the system using the ACUSON Freestyle mobile link app on a Microsoft® Windows® tablet*.

*Operates on a Microsoft® tablet device running Windows® operating system 8.1 or higher. For informational purposes only. Not intended for diagnostic use.

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

Standalone clinical images may have been cropped to better visualize pathology.

ACUSON Freestyle and Pixelformer are trademarks of Siemens Medical Solutions, Inc. USA.

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Siemens AG is under license.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare
Henkestraße 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
www.siemens.com/healthcare

Legal Manufacturer

Siemens Medical Solutions USA, Inc.
Ultrasound
685 E. Middlefield Road
Mountain View, CA 94043
USA
Phone: +1-888-826-9702

SIEMENS



Innovative to the Touch

Discover ultrasound insight powered by connectivity.

The ACUSON Freestyle Mobile Link App

Pairing the ACUSON Freestyle™ ultrasound system with Microsoft® Windows® devices, the ACUSON Freestyle™ mobile link app bridges the gap between the ultrasound insight you want and the equipment you use.

For more information, visit
[siemens.com/freestyle](https://www.siemens.com/freestyle)

Answers for life.



Simpler Workflow

The ACUSON Freestyle mobile link app creates a unique connection between you and your ACUSON Freestyle ultrasound system. By delivering valuable information directly to your Windows device, you have the freedom to work on the go.

Train smarter and farther than before

With the ability to securely share images and clips with your peers, the ACUSON Freestyle mobile link app supports your training and review of procedural techniques within and beyond your hospital's walls.

Get comprehensive ultrasound support—now

Easily share questions, concerns, or service requests in your Siemens query with the ACUSON Freestyle mobile link app's e-mail feature that automatically populates reference data—serial number, software version, IP address, etc.—from your ACUSON Freestyle ultrasound system.

The ACUSON Freestyle mobile link app is intuitively designed to optimize your processes by linking you directly to the powerful information on your ACUSON Freestyle ultrasound system. It's accessible innovation built specifically for you.

The products/features mentioned in this document may not be commercially available in all countries. Due to regulatory reasons their future availability cannot be guaranteed. Please contact your local Siemens organization for further details.

ACUSON Freestyle mobile link app operates on Microsoft® tablets running Windows® operating systems 8.1 or higher. For informational purposes only. Not intended for diagnostic use.




Effortless Understanding

The ACUSON Freestyle mobile link app capitalizes on some of today's most widely used mobile technologies, adding a sense of comfort and familiarity to your experience.

The ACUSON Freestyle mobile link app offers you a wide variety of accessible features, including:

- Easy, fast, and secure wireless setup
- Intuitive touch-screen interface
- Effortless controls for exploring your images and clips: just tap to select, swipe to navigate, stretch/pinch to enlarge/shrink



Seamless Connection

By delivering the information you want wherever you want it, the ACUSON Freestyle mobile link app assumes the role of digital ultrasound assistant. This integrated app streamlines your workflow by delivering powerful ultrasound insight that helps save time and allows for easier sharing of information.

Promote faster patient throughput

- The ACUSON Freestyle mobile link app offers in-app patient registration, allowing you to send your patient's information directly to your ACUSON Freestyle ultrasound system
- Simple, rapid, and secure access to your ultrasound information allows you to effortlessly attach images and clips to patient records—all from within the app

ACUSON Freestyle is a trademark of Siemens Medical Solutions USA, Inc.

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries.

Siemens Healthcare Headquarters
Siemens Healthcare GmbH
Henkestraße 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Legal Manufacturer
Siemens Medical Solutions USA, Inc.
Ultrasound
685 East Middlefield Road
Mountain View, CA 94043
USA
Phone: +1 888 826 9702
siemens.com/ultrasound